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C-Prime S plus

2. Composition

Modified acrylic acids, polyacrylic acid, methacrylates, catalysts, ethanol, stabilisators

3. Indication

C-Prime S plus is a one bottle light-cure single component priming and Dentin-Enamel adhesive, designed for bonding of composite / compomeres to natural Enamel and Dentin.

4. Contraindications

Placement is contraindicated if the patient is known to be allergic to any of the ingredients in C-Prime S plus.

5. Interactions

C-Prime S plus may not to be used together with underfillings that contains eugenol because phenolic substances could inhibit the polymerization process.

6. Warning

C-Prime S plus is sensitive to blue light, therefore the polymerisation can be started by daylight or operating light. For that reason light curing materials should be removed from the vials only immediately before application. Please note, that C-Prime S plus will not selfcure.

7. Application

- Isolation

Rubberdam is the recommended method of isolation.

- Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5 - 1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength. Cavity floor of deep excavations should be covered with a thin layer of hard setting Calciumhydroxide material.

- Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply C-Cid onto the enamel surfaces. Leave the etching gel in place for 30 - 45 seconds. After that apply the C-Cid additional onto the dentin surface and leave it for another 10 – 15 seconds. Rinse for 20 seconds with water. Dry it carefully in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of C-Prime S plus. The etched enamel bevel should have a chalky white appearance.

Etching precaution: It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. For the indications of compomeres etching is mostly not recommended, but it improves considerably the bonding with tooth material.

- Placement

Apply C-Prime S plus generously with a brush onto the enamel and dentin surfaces for at least 30 seconds with agitation. The dentin surface have to be wet in this time. Dry cautiously with oil free air for 15 seconds. Do not desiccate the dentin. Cure the C-Prime S plus coating by exposing its entire area with a dental halogen light unit ore strong LED light (e.a. Megalux) for at least 20 seconds .

- Placement of a second Layer of C-Prime S plus

Apply again the C-Prime S plus generously with a brush onto the adhesive surfaces for at least 30 seconds with agitation. Thoroughly dry again the primed areas with oil free air for 15 seconds and cure the C-Prime S plus coating by exposing its entire area to a dental halogen light unit for at least 20 seconds before placement of a light curing composite.

- Restorative Placement, Curing and Finishing

Refer to manufacturers instructions for placement, curing and finishing of the restorative material. After applying the composite on the wet surface C-Prime S plus develops together with the composite the optimal bonding values under oxygen-free conditions.

8. Storage stability

2 years

9. Storage and keeping instructions

Not above 25°C, refrigerated storage (+5°C to +8°C) is recommended.

10. Deliver units

5 ml C-Prime S plus ,
part of system packages

11. Date of issue

September 2009

12. Manufacturer and sales

MEGADENTA Dentalprodukte GmbH, D-01454 Radeberg, Germany
Tel. +49-3528-453-0 Fax +49-3528-453-21 Mail : info@megadenta.de