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**Directions of use**

# **C-Seal**

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## **Directions of use**

### **1. Information**

C-Seal light curing pit and fissure sealant

### **2. Ingredients**

C-Seal contains 50% (by weight) mineral and 0,1% Fluorine as Amino fluoride.  
The Methacrylate matrix contains Bis GMA, Bis EMA and BUDDMA.

### **3. Indication**

For sealing fissures and occlusal surfaces for the purpose of caries prophylaxis

### **4. Side effects**

The unpolymerized material may cause skin irritation.

Do not use this material on patients with a hypersensitivity against methacrylates.

### **5. Interactions**

C-Seal may not be used together with underfillings containing eugenol because phenolic substances could inhibit the polymerisation process.

### **6. Warning**

Avoid skin- and eyes contact with unpolymerized material. Since C-Seal cures with light, it is sensitive to blue light (operating light, sun light). For that reason bright light during application of C-Seal should be avoided.

### **7. Application**

-Cleaning and drying the teeth:

Before etching the tooth surface have to be thoroughly cleaned with an oil free polishing paste or a pumice powder. The surfaces must be dried.

-Enamel etching:

The cleaned enamel is etched with C-Cid. It is applied with the supplied brush. After 20-40s C-Cid is thoroughly rinsed with water and then carefully dried with oil free air. The etched and dried dentin enamel must show a smooth, chalky surface.

The etched enamel should not be contaminated by saliva or blood since this would lead to a deposit of protein that would inactivate the etched enamel surface. Should the etched enamel nevertheless become contaminated, then the etching process should be repeated.

-Application:

Apply C-Seal directly with the needle into the fissure. Application is recommended from the margin into the fissure, which allow to penetrate for 15-20 s. remove any air bubbles with the needle ore with a suitable instrument.

Light polymerisation with a blue light source ( e.a. Megalux) for 20 – 40 s. The surface of C-Seal has been carried out after polymerisation. The sealant parts of the teeth are easy to control with the white shading of C-Seal

### **8. Shelf life**

2 years

### **9. Precaution and Storage**

Refrigerated storage (+5°C to +8°C) is recommended; do not exceed 25°C.

**10. Delivery units**

single syringe C-Seal á 1.8 g  
assortment 4 x C-Seal á 1.8 g  
accessories.

**11. Date information prepared**

02/2008

**12. Manufacture and Sales**

MEGADENTA Dentalprodukte GmbH  
Carl-Eschebach-Strasse 1 A  
D-01454 Radeberg  
Germany

---

1340 0208 516