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## **1. Information**

Mega-A

Cold curing resin for repairs and relining of total and partial dentures

Shades: pink, blue-pink veined and clear

## **2. Composition**

Powder/liquid cold cure denture resin

Main component powder: polymethyl methacrylate

Component liquid: methyl methacrylate, dimethacrylate

## **3. Areas of application**

Resin for functional reconstruction, repairing, extending, supplementing as well as partial and total re-lining of dentures

## **4. Contraindications**

Skin contact should be avoided in patients with known hypersensitivity to methacrylates.

## **5. Side effects**

In individual cases, contact allergy may occur.

## **6. Interactions**

Unknown.

## **7. Hazard Information**

Read Safety data sheet first. When working with methyl methacrylates, direct skin and eye contact should be avoided. (Use personal protective equipment: protective gloves / eye protection). If swallowed, seek medical advice immediately. Ensure adequate ventilation during processing. Do not inhale vapors.

Mega-A liquid is highly flammable, the vapors may form explosive mixtures with air. Remove all sources of ignition. Close bottle immediately after usage.

## **8. Handling**

- Preparations

For embedding the prosthesis set up in wax, a semi-hard plaster is sufficient. After the wax has been scalded, the plaster is coated with MEGA-I or another alginate insulating agent.

*- Dosage and Handling*

Powder and liquid are intensively mixed in a ratio of **5 g powder : 3 ml liquid** within 20 s. The liquid is placed in the mixing vessel and the powder is interspersed until a low liquid level remains. The droppable liquid dough can be processed immediately. Apply MEGA-A evenly. After reaching the plastic consistency, do not touch and model again.

*- Polymerization*

The Curing time at 21 °C is approximately 13 min. Optimum polymerization takes place in the pressure pot at 45 °C at 2 bar (10 min). If processing without pressure curing unit, apply large amounts of MEGA-A in layers. Each layer must be completely cured and chilled before applying new dough.

*-Preparation and polishing*

After polymerization, the denture is removed from the model and processed with cross-cut milling cutters. Before final polishing with pumice stone, the denture is smoothed with abrasive paper of decreasing grain size. After production, MEGA-A contains an admissible residual amount of monomer < 2 % by weight. In order to prevent the risk of incompatibilities further, it is recommended to store the prosthesis in lukewarm water for at least 12 hours before it is inserted.

**9. Shelf life**

Liquid: 2 years

Powder: 3 years

**10. Precaution and storage**

Not above 25 °C

**11. Delivery units**

Powder: 80 g / 500 g / 1000 g

Liquid: 80 ml / 500 ml / 1000 ml

**12. Waste treatment**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

**13. Date information prepared**

November 2017

**14. Manufacturer and Sales**

MEGADENTA Dentalprodukte GmbH

Carl-Eschebach-Str. 1A / D-01454 Radeberg, Germany