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1. Information

Mega-F

Cold curing resin for denture bases and for model casting technology

Shades: pink, blue-pink veined and clear

2. Composition

Powder/liquid cold cure denture resin

Main component powder: polymethyl methacrylate

Component liquid: methyl methacrylate, dimethylacrylate

3. Areas of application

Resin for functional reconstruction of dentures, e.g. repairing, expanding, supplementing and partial and total relinings as well as for the preparation of dentures in the pressing process and in the model cast technique.

4. Contraindications

Avoid skin contact in patients with known hypersensitivity to methacrylates.

5. Side effects

In individual cases, contact allergy may occur.

6. Interactions

Unknown.

7. Hazard Information

Read Safety data sheet first. When working with methyl methacrylates, direct skin and eye contact should be avoided. (Use personal protective equipment: protective gloves / eye protection). If swallowed, seek medical advice immediately. Ensure adequate ventilation during processing. Do not inhale vapors.

Mega-F liquid is highly flammable, the vapors may form explosive mixtures with air. Remove all sources of ignition. Close bottle immediately after usage.

8. Handling

- Preparations

For embedding the prosthesis set up in wax, a semi-hard plaster is sufficient. After the wax has been scalded, the plaster is coated with MEGA-I or another alginate insulating agent.

- Dosage and Handling

Powder and liquid is mixed in a ratio of **10 g powder : 7 ml liquid**. The liquid is placed in the mixing jar and then the powder is quickly added and stirred immediately. A further addition of powder should only be carried out in small quantities and also quickly. Avoid the formation of blisters. Bubbles can be removed by gently swivelling. After mixing, MEGA-F can be poured for 2-3 minutes, after approx. 4 minutes, MEGA-F reaches the plastic processing phase of 3 to 3.5 minutes, during this period the dough can be modeled.

- Polymerization

After the total application time of 7 to 8 minutes the polymerization is carried out in the pressure vessel at a water temperature of 40-45 °C and 2 MPa (= 20 bar) within 6-7 min.

-Preparation and polishing

After polymerization, the denture is removed from the model and processed with cross-cut milling cutters. Before final polishing with pumice stone, the denture is smoothed with abrasive paper of decreasing grain size. After production, MEGA-F contains an admissible residual amount of monomer < 2 % by weight. In order to prevent the risk of incompatibilities further, it is recommended to store the prosthesis in lukewarm water for at least 12 hours before it is inserted.

9. Shelf life

Liquid: 2 years

Powder: 3 years

10. Precaution and storage

Not above 25 °C

11. Delivery units

Powder: 80 g / 500 g / 1000 g

Liquid: 80 ml / 500 ml / 1000 ml

12. Waste treatment

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

13. Date information prepared

November 2017

14. Manufacturer and Sales

MEGADENTA Dentalprodukte GmbH

Carl-Eschebach-Str. 1A / D-01454 Radeberg, Germany