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## 1. Information

Mega-F

Cold curing resin for dentures and for precision-cast technique

Shades : pink, blue-pink veined and clear

## 2. Composition

Powder/liquid cold cure denture resin

Main component powder: Methylmethacrylate, Copolymer

Component liquid: Methylmethacrylate

## 3. Indications

Dental denture repair resin for functional reconstruction, repair, partial and total re-lining of dentures

## 4. Contraindication

Allergic reactions in wearers to dental prosthesis base resin are rare. Placement of Mega-F is contraindicated if the patient is known to be allergic to any of the ingredients in Mega-F.

## 5. Side Effects

none admits

## 6. Warning

Mega-A liquid contains high flammable Methylmethacrylate. Irritating to eyes, respiratory system and skin. Avoid contact with skin, wear suitable gloves. Keep away from sources of ignition ( no smoking) , keep out of reach of children. If swallowed, seek medical advice immediately and show the product label.

## 7. Application

### - Preparations

Use hard plaster for model. Remove wax residue thoroughly from the resin parts. Insulate watered plaster model twice with a thin layer of MEGA-I separating solution (ore similar).

### - Dosing and Handling

Powder and liquid is mixed at a ratio of 10 g and 7 ml. The liquid is placed into the mixing jar and the approximate required quantity of powder is added quickly and stirred immediately. Additional powder should be added only in small quantities and also quickly. Avoid air entrapping, swivel mixing jar to remove same. The MEGA-F dough is 2 -3 min pourable after mixing. MEGA-F has a plastic handling phase after approx. 4 min which lasts for about 3 to 3.5 min. During this period the dough is stable and can be moulded.

### - Polymerization

Curing in a pressure polymerization unit starts between 5th min and 8 min, the curing time in a pressure unit is 6 – 7 min, the water temperature 45°C and the pressure 2 bars.

## 8. Shelf life

Liquid: 2 years

Powder: 3 years

## 9. Precautions and Storage

Not above 25°C

## 10. Deliver units

Powder: 80 g / 500 g / 1000 g

Liquid : 80 ml / 500 ml / 1000 ml

## 11. Date information prepared

September 2010

## 12. Manufacture and sales

MEGADENTA Dentalprodukte GmbH

D-01454 Radeberg, Germany