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Direction of use

Mega- L

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Directions of use

1. Information

Mega-L
Heat-curing resin for dentures in form of powder / liquid

2. Composition

Powder/liquid hot-cure denture resin
Main component powder: Methylmethacrylate, Copolymer
Component liquid: Methylmethacrylate

3. Indications

For dental dentures manufactured by the press process.

4. Contraindication

Allergic reactions in wearers to dental prosthesis base resin are rare.
Placement of Mega-L is contraindicated if the patient is known to be allergic to any of the ingredients in Mega-L.

5. Side Effects

none admits

6. Warning

Mega-L liquid contains high flammable Methylmethacrylate. Irritating to eyes, respiratory system and skin. Avoid contact with skin, wear suitable gloves. Keep away from sources of ignition (no smoking) , keep out of reach of children. If swallowed, seek medical advice immediately and show the product label.

7. Application

Recommended mixing ratio is 3 g powder : 1 ml liquid. Powder and liquid can be freely dosed. First place the liquid in the mixing vessel, and then quickly shake in powder until a thin, dry layer of powder remains lying on the surface. Tip off this powder surplus after 30 s. Then thoroughly mix with a broad spatula. The relative thin consistency of the dough should not lead one to add further powder. Allow the dough to swell in the mixing vessel for about 5-10 min. At normal room temperature, MEGA-L will then remain at processing consistency for up to 60 min,

Filling and pressing the flask

Slightly overfill the flask with dough. Gradually close the flask under a press in intervals, and leave for about 10 min under final pressure.

In the case of thick work pieces add small supplementary quantities of dough.

The temperature of the flask when it is packed should be between roomtemperature and 45 °C

Polymerization

MEGA-L can be polymerized by the short-time, normal or long-time method.

Water bath polymerization

Place the flask in at least 3 l of luke-warm water. Gradually heat the water bath to 70°C, and maintain this temperature in accordance of the thickness for approximately 30 – 60 min. Than bring the water bath to the boil. Allow to boil for 30 min.

Short-time polymerization in a water bath

A prerequisite to short-time polymerization is the maintenance of an internal flask temperature of 80°C for at least 15 min. The flask should therefore be placed in boiling water. Wait until the water returns to the

boil, and then switch off the heat immediately and cover the bath. The flask must remain in the water bath for at least 45 min, and then be allowed to cool on the air. At least 2 l of water per flask should be used. If only 1 flask is polymerized, at least 3 l are necessary.

Important:

The properties of MEGA-L are only then ensured in short-time polymerization if the following points are noted:

The above water quantities must be used to ensure a sufficient heat reserve. The water bath should not be left to boil after returning to the boiling point, particularly when polymerizing thick workpieces. This avoids pores and bubbles which could occur as a result of overheating. If short-time polymerization is carried out using an electric hotplate, remove the water bath from the hotplate. The flask must remain in the water bath for at least 45 min.

Long-time polymerization in a water bath

Place the flask in a polymerization unit which has been filled with cold water. One then polymerizes at 80°C for 10 hours in this unit. The long polymerization period ensures that the temperature rises gradually.

Polymerization in a hot-air cabinet

The material is polymerized at a temperature of 90°C for 10 hours. This ensures that the polymerization temperature rises gradually. Slow cooling of the flask is also necessary in this polymerization method.

8. Shelf life

Liquid: 2 years

Powder: 3 years

9. Precautions and Storage

Not above 25°C

10. Deliver units

Powder: 80 g / 500 g / 1000 g

Liquid : 80 ml / 500 ml / 1000 ml

11. Date information prepared

January 2003

12. Manufacture and sales

MEGADENTA Dentalprodukte GmbH
D-01454 Radeberg, Germany