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Directions of use

N-Seal

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Directions of use

1. Information

N-Seal light curing pit and fissure sealant with Nano-fillers

2. Ingredients

N-Seal contains 60% (by weight) mineral and 0,1% Fluorine as Amino fluoride. The Methacrylate matrix contains Bis GMA, Bis EMA and BUDDMA.

3. Indication

For sealing fissures and occlusal surfaces for the purpose of caries prophylaxis

4. Side effects

The unpolymerized material may cause skin irritation.

Do not use this material on patients with a hypersensitivity against methacrylates.

5. Interactions

N-Seal may not be used together with underfillings containing eugenol because phenolic substances could inhibit the polymerisation process.

6. Warning

Avoid skin- and eyes contact with unpolymerized material. Since N-Seal cures with light, it is sensitive to blue light (operating light, sun light). For that reason bright light during application of N-Seal should be avoided.

7. Application

-Cleaning and drying the teeth:

Before etching the tooth surface have to be thoroughly cleaned with an oil free polishing paste or a pumice powder. The surfaces must be dried.

-Enamel etching:

The cleaned enamel is etched with C-Cid. It is applied with the supplied brush. After 20-40s C-Cid is thoroughly rinsed with water and then carefully dried with oil free air. The etched and dried dentin enamel must show a smooth, chalky surface.

The etched enamel should not be contaminated by saliva or blood since this would lead to a deposit of protein that would inactivate the etched enamel surface. Should the etched enamel nevertheless become contaminated, then the etching process should be repeated.

-Application:

Apply N-Seal directly with the needle into the fissure. Application is recommended from the margin into the fissure, which allow to penetrate for 15-20 s. remove any air bubbles with the needle ore with a suitable instrument.

Light polymerisation with a blue light source (e.a. Megalux) for 20 – 40 s. The surface of N-Seal has been carried out after polymerisation. The sealant parts of the teeth are easy to control with the white shading of N-Seal

8. Shelf life

2 years

9. Precaution and Storage

Refrigerated storage (+5°C to +8°C) is recommended; do not exceed 25°C.

10. Delivery units

single syringe N-Seal á 1.8 g
assortment 4 x N-Seal á 1.8 g
 accessories.

11. Date information prepared

02/2008

12. Manufacture and Sales

MEGADENTA Dentalprodukte GmbH
Carl-Eschebach-Strasse 1 A
D-01454 Radeberg
Germany

1340 0208 517