# T-Fill eco

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## Direction for use

### 1. Product name

T-Fill eco Temporary Filling material of zinc oxide base

### 2. Composition:

powder: zinc oxide, pharmaceutically quality,

liquid: water-solution with zinc sulphate and acid solution

### 3. Indication

T-Fill eco is used for temporary fillings and for temporary fixation of crown- and bridges.

### 4. Contraindications

So far unknown

### 5. Interactions

None reported

### 6. Application:

In all classes of cavities this may be kept to the minimum required for caries removal. Undercuts and dovetails are possible for fixation but not required. After the cavity is cut it has to be washed with water and dried with a gentle stream of air a desiccation is to avoid. Deep cavities or possible pulp exposure areas should be lined with a calcium hydroxide material.

The mixing ratio of T-Fill eco powder/liquid 1:1. The minimum that should be used is scarcely 1 scoops of powder and 1 drops of liquid. The mixing should be done on the mixing pad or a clean glass slab. The powder is loosen up by stirring or shaking the bottle several times, then placed onto mixing pad and divided into equal portions. The required number of liquid drops is dispensed onto the mixing pad. With a plastic or stainless steel spatula is one of the three portions of powder incorporated onto the liquid. In successive portions will be mixed until T-Fill eco reaches a homogenous creamy consistency. The mixing should be completed in 30 s. Rinse spatula and slab immediately after use with water.

#### Mixing and Setting Times(23°C)

Mixing rate powder/liquid: 2:1

Mixing Time: 30 s

Working Time: at least 1 min

Initial Setting Time: 3 min

Total setting time 5 min

T-Fill eco can be removed from cavity with instruments ore cutting bores.

### 7. Storage stability

2 years

### 8. Storage and keeping instructions

dry, not above 25°C

### 9. Administration and package size

bottle of 80 g powder / 80 ml liquid

accessories

### 10. Date information published

04/2011

### 11. Manufactures and sales

MEGADENTA Dentalprodukte GmbH

D-01454 Radeberg, Germany

Internet: [www.megadenta.de](http://www.megadenta.de) / [info@megadenta.de](mailto:info@megadenta.de)